

**NEW!**

**FDA Cleared TDM for Infliximab & Adalimumab**

Take control of TDM today with results in less than 5 minutes

ProciseDx™ provides immediate serum drug level results of Adalimumab and Infliximab in 5 minutes, right from your laboratory. Give your Crohn's and Colitis patients same day TDM results so you can take control of patient management today, not days or weeks later



- ✓ FDA Cleared TDM for Infliximab and Adalimumab
- ✓ Easy to perform with TDM results in less than 5 minutes
- ✓ Get important drug concentration information without compromising lab-quality results
- ✓ Fits any moderate complexity clinical lab workflow

**Ask for a demo at ProciseDx Booth 330**

**Product Information**

|                         | Procise ADL ✓                      | Procise IFX ✓                                     |
|-------------------------|------------------------------------|---|
| <b>Sample Type</b>      | <b>Serum</b>                       | <b>Serum</b>                                      |
| <b>Reportable Range</b> | 1.0 µg/mL - 50.0 µg/mL             | 1.2 µg/mL - 50.0 µg/mL                            |
| <b>Time to Results</b>  | <b>Less than 5 minutes</b>         | <b>Less than 5 minutes</b>                        |
| <b>Stability</b>        | 2 years at room temperature        | 2 years at room temperature                       |
| <b>Validated for</b>    | Adalimumab (Humira® and Amgevita®) | Infliximab (Remicade®, Inflectra® and Renflexis®) |

Calibrated against WHO Standard